

BLOOD SAMPLE COLLECTION INSTRUCTIONS

Do NOT store kits in an area where the temperature range is outside 18–30°C or 65–86°F.


Do NOT expose blood to temperatures outside the range of 43–98°F or 6–37°C.

1. Complete requisition form

- Select the appropriate Signatera screening option.
- Complete ALL appropriate fields on the requisition form.

2. Collect blood sample


One 6mL EDTA tube
AND two 10mL DNA
Streck tubes

- Fill two Streck tubes for all patients. An additional 6 mL EDTA is required for first time patients only.
- Fill all tubes completely. If insufficient volume is obtained, draw an additional tube.
- Allow 60–90 seconds for each tube to fill.
- Use 21 gauge straight needle.
- Vein collapse may require second venipuncture with a fresh tube.
- **Gently** invert tube at least 8 times immediately after draw in order to thoroughly mix blood with reagents.
- Do **NOT** shake vigorously.

3. Pre-pack sample


A
Complete information
on tube labels.

B
Place filled and labeled
tubes into absorbant
sleeve and into the metallic
envelope. Place the
metallic envelope into the
biohazard bag.

4. Shipping

A
Pack the blood sample, the completed requisition form and the cool pack included, into the Signatera Blood Collection Kit box.
Return* the box to the laboratory of Unilabs Genetics in Lausanne at the address indicated below.

Contact details of lab/courier service:
Dr Kostas Nikopoulos/Dr Mattia Schmid
Diagnostic moléculaire
Rue de la Vigie 5 - CH - 1003 Lausanne
T: 021 321 40 51 ou 52 ou 55
Samples are accepted:
Monday-Thursday: 09h00-16h00*

*We highly encourage you to contact us prior to sampling in order to further personalise and optimise the sampling conditions.

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Natera, Inc.
Qarad B.V.